


JOBST® FarrowWrap® 4000 A-D Below Knee

Ready-to-wear order form and prescription request

Email order to Customer Services on: **compression.uk@jobst.com**
 Fax order to Customer Services on: **0345 122 3450**

Date: _____ Purchase Order No.: _____ Patient Name: _____ DoB: _____
 Measured By: _____ Tel: _____ Email: _____
 Delivery Address: _____ Invoice Address: _____

Legpiece



Measuring instructions

Legpiece

1. Measure narrowest ankle **B** and widest calf **C** circumferences
2. Measure length posteriorly (back of leg) following the contour from **B** (narrowest part of the ankle) to **D** (two finger widths below the crease of the knee)
3. Determine size (**XS-L**) & length (**Regular or Tall**) using chart below

When ordering JOBST FarrowWrap 4000, 1 x JOBST FarrowHybrid (Liners) wide is also provided. See table below for the JOBST FarrowHybrid size that will be provided with each JOBST FarrowWrap 4000 size ordered. If additional or alternative JOBST FarrowHybrid (Liners) sizes are required, please order separately

JOBST FarrowHybrid (Liners)

1. Measure **A1** circumference
2. Record **shoe size**
3. Determine size (**S-XL**), width (**Regular or Wide**) & colour (**Tan or Black**)

Check out our measuring guide videos

Scan the QR code for help with measuring for JOBST FarrowWrap



JOBST® FarrowWrap® 4000 Legpiece

Size: X Small Small Medium Large

Length **B-D** Regular (30-33 cm) Tall (34-37 cm)

Colour: Tan Black

PIP Code: -

Product Code: -

Quantity Required:

Repeat prescription required every months

JOBST® FarrowHybrid® (Liners)

Size: Small (2-5) Medium (6-8) Large (9-11) X Large (12-13) (shoe size)

Width **A1** Regular (20-27cm) Wide (28-35cm)

Colour: Tan Black (available in Wide variant only)

PIP Code: -

Product Code: -

Quantity Required:

Repeat prescription required every months

	X Small	Small	Medium	Large
JOBST FarrowWrap 4000 Legpiece				
C Calf	25-38cm	32-48cm	37-56cm	43-67cm
B Ankle	17-23cm	19-28cm	25-37cm	29-43cm
B-D Regular	30-33cm	30-33cm	30-33cm	30-33cm
B-D Tall	34-37cm	34-37cm	34-37cm	34-37cm

NOTE: for **B-D** length measure posteriorly (back of leg) and follow contours

	JOBST FarrowWrap 4000 size ordered			
	X Small	Small	Medium	Large
JOBST FarrowHybrid size received with wrap				
Regular length	Medium-wide	Medium-wide	Medium-wide	Medium-wide*
Tall length	Large-wide	Large-wide	Large-wide	Large-wide*

*Do not apply JOBST FarrowHybrid if calf circumference exceeds 60cm. When a large JOBST FarrowWrap 4000 is ordered, an additional non-compressive liner will be provided to accommodate calf circumference greater than 60cm. Consider a JOBST FarrowWrap footpiece in this circumstance.

JOBST FarrowWrap 4000 Order Codes

JOBST® FarrowWrap® 4000 Legpiece (30-40mmHg)

For patients with **mild to moderate** oedema, with or without a leg ulcer

Product Code	Size	Length	Colour	PIP Code
76661-00	Extra small	Regular	Beige	408-1543
76661-08	Extra small	Regular	Black	408-1584
76661-04	Extra small	Tall	Beige	408-1501
76661-12	Extra small	Tall	Black	408-1469
76661-01	Small	Regular	Beige	408-1550
76661-09	Small	Regular	Black	408-1592
76661-05	Small	Tall	Beige	408-1519
76661-13	Small	Tall	Black	408-1477
76661-02	Medium	Regular	Beige	408-1568
76661-10	Medium	Regular	Black	408-1600
76661-06	Medium	Tall	Beige	408-1527
76661-14	Medium	Tall	Black	408-1485
76661-03	Large	Regular	Beige	408-1576
76661-11	Large	Regular	Black	408-1618
76661-07	Large	Tall	Beige	408-1535
76661-15	Large	Tall	Black	408-1493



JOBST® FarrowHybrid® (20-30mmHg)

For the treatment of **mild to moderate** oedema in the foot and ankle

Product Code	Size	Colour	PIP Code
76664-00	Small - Regular (shoe size 2-5)	Tan	408-1626
76664-01	Medium - Regular (shoe size 6-8)	Tan	408-1634
76664-02	Large - Regular (shoe size 9-11)	Tan	408-1642
76664-03	Extra Large - Regular (shoe size 12-13)	Tan	408-1659
76664-12	Small - Wide (shoe size 2-5)	Tan	408-1667
76664-13	Medium - Wide (shoe size 6-8)	Tan	408-1675
76664-14	Large - Wide (shoe size 9-11)	Tan	408-2384
76664-15	Extra Large - Wide (shoe size 12-13)	Tan	408-1683
76664-16	Small - Wide (shoe size 2-5)	Black	408-6765
76664-17	Medium - Wide (shoe size 6-8)	Black	408-6757
76664-18	Large - Wide (shoe size 9-11)	Black	408-6740
76664-19	Extra Large - Wide (shoe size 12-13)	Black	408-6732



Do not apply JOBST FarrowHybrid if calf circumference exceeds 60cm, consider a JOBST FarrowWrap footpiece in this circumstance

Find ordering and prescribing information fast with JOBST Code Finder

No more searching through a catalogue, use our web app on your mobile device to find the information you need.

Scan the QR code or search for [codefinder.jobst.com](#) on your mobile device.



Fast and free home delivery of prescription garments with JOBST Delivered*

To use this dispensing service, send this form and prescription to: **JOBST Delivered, FREEPOST, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG** or order online at [order.jobst.com/UK](#) and select JOBST Delivered as the dispenser (note: prescription to be posted to JOBST Delivered unless using EPS).

If the patient is using an Electronic Prescription Service (EPS), select code **FJ708** to nominate Daylong Direct as their dispenser of choice.

For envelopes, please contact: [jobstdelivered@daylongdirect.co.uk](#) or **0808 196 6250**

*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG
Privacy Policy: [www.daylong.co.uk](#)

Go digital and order using JOBST Online

Bringing you an easier way to order JOBST. Order direct with payment by account or on prescription by selecting JOBST Delivered* as the dispenser on JOBST Online.

Scan QR code to register or visit [order.jobst.com/UK](#)



*JOBST Delivered is operated by Daylong Direct, 10 Cossall Industrial Estate, Ilkeston, Derby, DE7 5UG

Privacy Policy: [www.daylong.co.uk](#)